

Online practical gynecological workshops October 25, 2023

Urological tapes are knitted of transparent monofilament polypropylene yarn. On both endings the tape is equipped with blue handles which make it possible to fix it on an appropriate applicator. The handles are also made of monofilament polypropylene yarn. Each handle has additional protection against pulling the tape out with the use of force necessary for putting it through anatomical structures. The tape characterises with: biostability, thanks to which it does not lose physical properties while long-term contact with tissues and systemic fluids, monofilament structure which minimises the risk of bacterial infection, and optimally elaborated spatial structure, where the size of pores makes it possible for quick overgrowth with connective tissue.

Application:

- for surgical treatment of stress urinary incontinence in women with the use of tension-free TVT or TOT method


Types:

- sterile (sterilised with ethylene oxide)

Size and packing:

- length of 45 cm and 60 cm
- packed in double paper-film packaging with gusset and then in a cardboard box by 1 pc

More information on: www.dallop-nm.com



MD

Sterile devices are certified by the Notified Body PCBC SA and bear the CE mark with number 1434

>>