



MATOSSET PICTOGRAMS - INSTRUCTION

How to read out the information - read horizontally:



Description of particular product parameters:

Pictogram	Explanation	Pictogram	Explanation
Quantity		Product type	
1 (2 x 10)	1 piece tied: 2 piles by 10 pieces		drape: of nonwoven, with adhesive, with opening, with opening and adhesive, under the buttocks
Product type			T-shaped with opening, for ophthalmology, for angiography
	swabs, swabs with Y and O cut		for urology, with „U”shape cut-out
	neurosurgical swab		for laparoscopy, for craniotomy, for cystoscopy
	tampons: ball, bean, whisker, triangle, thimble, grain, pepper,		vertical drape, Mayo table cover
	laryngological		instrument table cover, 5-chamber liquid pouch, Seni Soft absorbent underpad, film pouch for organ
	abdominal swab with cotton loop, ribbon		liquid pouch, medical fixing tape / hook and loop fixing tape
	Absorgyn gynaecological towel		
	Absora: absorbent swab, high-absorbent swab, tamponade dressing, cotton-wool and gauze tampon		
	Absorba: eye dressing, nose swab		
	Tubula tubular bandage, bandage with a clip, Matosoft orthopaedic padding		
	wound closure strips, adhesive tape with perforation, tourniquet		
	adhesive tape on a bobbin, Medisorb F dressing		

Pictogram	Explanation	Pictogram	Explanation
Product type		Product type	
	<p>hand towel, Fixopore S nonwoven adhesive dressing, Fixopore S eye dressing, film-nonwoven dressing for cannula fixation</p> <p>leg cover, arm cover, lamp handle cover</p> <p>round cover for medical equipment, rectangular cover for medical equipment, C-arm cover, Camera Cover for medical device cables</p> <p>surgical gown: Standard, Perfect, apron, sleeve protectors, tied shoe covers, cap</p> <p>nurses cap, surgical cap, surgical cap with sweatband, surgical cap with ties</p> <p>examination glove, surgical glove, face mask with ties, with elastic band</p> <p>container: rectangular, round, tray</p> <p>container: 2-chamber, 3-chamber, kidney dish</p> <p>container: 3-chamber, for cables 2500ml, safe box for scalpels, safe box for scalpels and needles</p> <p>sponge stick, cotton swab, tongue depressor</p> <p>eye spears for ophthalmology, plastic scissors sharp-blunt, Littauer scissors</p> <p>scissors: sharp-blunt, sharp-sharp, for episiotomy</p> <p>scissors, umbilical scissors, clamps</p> <p>sponge forceps, Kocher forceps/clamps, needle holder</p> <p>Halsted Mosquito forceps, Pean forceps straight, Backhouse clamp</p> <p>umbilical cord clamp cutter, skin staples remover</p>		<p>blade no. 10, no. 11, no. 15</p> <p>blade no. 20, no. 24</p> <p>stitch cutter, scalpel</p> <p>safe scalpels</p> <p>anatomical forceps, umbilical cord clamp, surgical forceps</p> <p>syringe, needle, drain</p> <p>hydrodelineation needle 0.5x22mm, olive cannula 0.6x22mm, curved irrigation cannula</p> <p>Foley catheter, urine bag, organ bag</p> <p>stitches, connector, 3-way stopcock with tubing</p> <p>transfusion apparatus, venflon, non-return valve</p> <p>Redon bottle, monopolar handle, Codan line</p> <p>Scalpfix scalp clip cartridge, 3-way stopcock</p> <p>barraquer, barraquer with closed loops, drape clamp, neurosurgical swab</p> <p>Introcan safety needle</p> <p>marker, baby measuring band, electrode cleaner</p> <p>3-way ramp of stopcocks</p>

Pictogram	Explanation	Pictogram	Explanation
Symbols			
<p>13,17, 20 </p> <p>30 g/m² , 40 g/m² </p> <p> </p> <p>4 </p> <p>SMS, BP, L2, L3, TF, FB, PE</p> <p>  powder-free powdered</p>	<p>thread count for gauze</p> <p>basis weight of nonwoven</p> <p>gauze/nonwoven with X-ray detectable thread</p> <p>number of plies fo gauze/ nonwoven products</p> <p>symbols of nonwoven/laminate/film types</p> <p>information on latex/powder presence/absence</p>		
Size			
<p>320x180 (Ø8)</p> <p>1 m x 2 cm</p> <p>200 ml</p> <p>nr 11</p> <p>CH14 70cm</p> <p>M, L, XL, XXL</p> <p>6; 6,5; 7; 7,5; 8; 8,5; 9</p> <p>XS, S, M, L, XL</p>	<p>320 cm x 180 cm, opening of 8 cm diameter (same units of measure)</p> <p>1 m x 2 cm (different units of measure)</p> <p>capacity 200 ml</p> <p>numbering used for describing blades</p> <p>sort and lenght of drain</p> <p>sizes of surgical gowns</p> <p>sizes of surgical gloves</p> <p>sizes of examination gloves</p>		
2/IJ-08-001-12.01.2023			