


SYMBOL	MEANING
	CE with number – confirmation of meeting EU requirements
	medical device
	manufacturer
	date of manufacture
	reference number of the product
	sterile product; sterilised using steam
	sterile product; sterilised using ethylene oxide
	sterile product; sterilised using irradiation
	sterile product; no information about the type of sterilisation when more than one method can be used
	non sterile product
	batch code – the manufacturer’s code
	use-by-date – indicates the date after which the medical device is not to be used
	do not re-use – product for single use only
	do not re-sterilise
	do not use if package is damaged
	single patient – multiple use
	not allowed for food contact
	allowed for food contact
	keep away from sunlight
	keep dry
	temperature limits to which the medical device can be safely exposed
	product contains natural rubber latex
	product is latex-free

	product gives protection against microorganisms
	product gives limited protection against chemicals
	instructions for use attached to the product
	range of humidity to which the medical device can be safely exposed
	„Green dot” means that the manufacturer of the product contributes to the cost of recovery and recycling
	information about the raw material the given packaging is made of
	throw away your trash
	non-pyrogenic product
	hand wash
	wash at the indicated temperature
	drip dry
	do not iron
	cannot be bleached with chlorine
	do not dry clean
	wash at normal temperature in the tumble drier
	must not be bleached
	not for contact with food
	indicator showing from which side should the packaging be opened
	this side up
	handle with care
	quality control approval
	registration number for Ukraine - informs since when the product is registered
	registration number for Russia - informs since when the product is registered

Per. № ФСЗ 2012/12667
от 22.08.2012 г.