

---

# Matoset Matodrape

QUALITY COMPREHENSIVITY SAFETY



---

## OPERATING BLOCK

---



 TZMO SA

---

## TZMO'S MEDICAL BRANDS

---



**TRICOMED**

Ladies and Gentlemen!

Medical devices under the **Matopat** brand bring almost 70 years of tradition and history, and above all, the development of a wide range of products intended for the operating theatre. This development is possible due to the constant introduction of cutting edge technologies, as well as close cooperation with medical specialists and authorities, and thanks to the flexibility of the offer.

All these activities are aimed at creating new products, gradual improvement of the current offer and meeting the expectations of our customers. Thanks to this, we design products that reflect the actual needs of medical staff.

When designing our products, we closely cooperate with the **Matopat Specialist Hospital**, where in the very centre of the hospital, in the operating theatre, we test our ideas and we use feedback from specialists regarding the parameters of our products. The result of this cooperation are **Matodrape** surgical drapes and **Matoset** operative sets, which allow for comprehensive preparation and carrying out the operation.

At the turn of the century, the model of functioning of health care facilities has changed. To be able to focus on their main mission - treating patients - hospitals increasingly decide to outsource processes not directly related to the treatment of patients. To meet these trends, **TZMO** has been building a network of **CitoNet Service Centres** for several years, which, using the washing and sterilisation processes, provide our customers with clean bedclothes, sterile surgical instruments and sterile operating clothing as well as sterile sets consisting of dressing material, drapes, gowns and surgical instruments. Today, the **CitoNet Service Centre** network consists of six locations in Poland: Łódź, Bydgoszcz, Kraków, Szczecin, Warsaw and Swarzędz. The management of the above centres are carried out in accordance with the highest standards, based on international **ISO 9001**, **ISO 14001** and **ISO 13485** standards, and the disinfection, sterilisation and washing processes are validated, thanks to which we have certainty that their effectiveness is repetitive.

A handwritten signature in blue ink, reading "Małgorzata Torusiel".

TZMO Director  
for Sales and CitoNet Services

## contents


---

	HEAD AND NECK SURGERY	7-10
	VASCULAR SURGERY	11-16
	GENERAL SURGERY	17-23
	GYNAECOLOGY AND OBSTETRICS	24-31
	OPHTHALMOLOGY	32-34
	ORTHOPAEDICS AND TRAUMATOLOGY	35-44
	UROLOGY	45-50
	CLOTHING	51-54

## legend of markings

---

### ■ Colours of components:

	----- surgical drapes
	----- covers and auxiliary elements
	----- clothing
	----- dressing materials
	----- biomaterials
	----- instruments, accessories

### ■ The raw material marking of the main drape can be found in the product code, e.g.:

MA-991-ZOTM-013	----- SMS nonwoven
MA-991-ZOBP-001	----- BP cotton-like nonwoven
MB-991-ZOTF-003	----- TF laminate
MB-991-ZOL2-042	----- L2 laminate
MB-136-ZOL3-010	----- L3 laminate
MB-135-ZOMI-008	----- made of various raw materials (MIX)
MA-136-ZOPE-001	----- PE film

### ■ Packaging types:

Blister pack (BLI)	----- soft blister packaging made of medical paper and transparent film
Paper-film pouch (PFP)	----- pouch made of medical paper and transparent film
Tyvek®-film (Tyvek®)	----- very durable packaging for packing heavier sets made of synthetic Tyvek® material and transparent film



■ surgical drapes

■ covers

■ elements organising work:

- fixing tapes
- organisers
- (liquid) pouches
- waste bags

■ surgical clothing:

- gowns
- sets
- gloves
- caps
- face masks

■ dressing materials

■ biomaterials

■ surgical instruments:

- made of metal and plastic

■ auxiliary elements:

- containers
- syringes and needles
- drains
- bottles





QUALITY  
**COMPREHENSIVITY**  
SAFETY

**Matoset** operative set

---

All the elements necessary to perform a given surgery in one complete set.

# Matodrape and Matoset sets are modern solutions to improve efficiency in the operating room.

The sets are:

- sterile, which minimises the risk of nosocomial infections,
- disposable, which eliminates the costs associated with repeated washing, ironing and sterilisation - elements typical for reusable products,
- custom-made - composition adapted to a given procedure, which reduces the burden on medical personnel with the need to complete individual elements, and thus shortens the time to prepare for the procedure,
- complete, which reduces the number of intermediate packages; this ensures greater safety and reduces the amount of waste.

## Matodrape

sets composed of surgical drapes,  
accessories and clothing



## Matoset


sets might contain:  
drapes, accessories and clothing  
as well as instruments and dressing materials



**Barrier\***  
laminates used  
in surgical  
drapes:



**TF laminate**  
2-layer, durable, absorbent  
and fully barrier



**L2 laminate = Blue Special**  
2-layer, very durable, very  
absorbent and fully barrier




**L3 laminate = Blue Comfort**  
3-layer, the most durable,  
highly absorbent and fully barrier

\* barrier to fluid penetration,  
microorganisms  
in accordance with EN 13795

plus

**PLUS** – an additional absorbent layer with exceptional absorption properties used in the critical area of the main drape perfectly distributes blood, liquids used during the procedure and body fluids. The fluids are distributed in all directions, even on a vertical surface. An additional absorbent layer allows to increase the total absorbency of the drape while leaving it soft and easy to drape in the surgical field.



## Handy TAG labels

**Matodrape**  
Zestaw do angiografii nr 3 CF

Do not use if package is damaged...  
LOT TEST 2018-05

MD

Manufacturer / Hersteller / Fabricante / Wytwórca:  
TIZMO S.A. ul. Żółkiewska 20/22, 01-100 Toruń, Polska

product group

set name

set composition

additional information

technological symbols

4 stickers (TAGs) for placing in medical records

producer

**Matoset**  
REF MA-991-ZOL2-661

Ophthalmology set  
Ophthalmologie-Set  
Trousse interventions ophtalmologiques  
Komplektas oftalmologinėms operacijoms  
Komplektas acu operacjăm  
Комплеккт для офтальмологических операций  
Zestaw do operacji okulistycznych

1 150x150  
1 150x120

1 17 2=5x6  
12 17 8=5x5  
2 LL 3ml  
1 LL 5ml  
1 LL 20ml  
2 60x40  
1 30ml

L / 125cm STANDARD  
M / 115cm STANDARD

Do not use if package is damaged...  
LOT TEST 2020-01

CE 2274 STERILE EO MD

Manufacturer / Hersteller / Fabricante / Wytwórca:  
TIZMO S.A. ul. Żółkiewska 20/22, 01-100 Toruń, Polska

## Practical double packing with a dispenser option



The catalogue includes our standard offer. The main drape of the set is usually available in a **TF, Blue Special and Blue Comfort** version. It is also possible to create a custom-made sets to meet individual needs with a wide range of components. Please contact our Sales Representative.

# CRANIOTOMY




plus



## Matoset

	MB-991-ZOL3-028	Blue Comfort	Tyvek®	pieces
A -	Drape 300cm×240cm with opening 30cm×20cm, incise film, liquid collection pouch, absorbent layer			1
B -	Drape 60cm×50cm with adhesive tape			4
C -	Instrument table cover 190cm×150cm			1
D -	Mayo table cover 145cm×80cm			1
E -	Adhesive Velcro tape 2cm×25cm			1
F -	Nonwoven hand towel 40cm×20cm			2
G -	Lamp handle cover 12cm×14cm			2
H -	Gauze swab 17-thread 16-ply with X-ray detectable thread 7.5cm×7.5cm			50
I -	Gauze swab 17-thread 8-ply 10cm×20cm			8
J -	Nonwoven swab 4-ply 40g/m <sup>2</sup> with X-ray detectable thread 7.5cm×7.5cm			50
K -	Neurosurgical swab of 100% cotton 6-ply with X-ray detectable thread 1.5cm×15cm – 10 pieces on a cardboard			1
L -	Gauze abdominal swab 17-thread 6-ply with X-ray detectable thread 45cm×45cm			1
M -	Sponge for cleaning electrodes			1
N -	Plastic container 700ml			1
O -	Scalpel blade disposal box with adhesive tape			1
P -	Syringe 2-piece L 20ml			2
Q -	Scalpel blade no. 11			1
R -	Scalpel blade no. 24			1
S -	Urine bag 2l			1
T -	Scalp skin suturing clips Scalpfix (cartridge by 10 pcs)			2




# CRANIOTOMY



## Matodrape

<b>MB-136-ZOL2-005</b>	Blue Special	PFP	pieces
------------------------	--------------	-----	--------

A -	Drape 240cm×200cm with opening 20cm×25cm, <b>incise film</b> and liquid collection pouch*	1
B -	Drape 60cm×50cm with adhesive tape	2
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Medical fixing tape 9cm×50cm	2
F -	Nonwoven hand towel 40cm×20cm	2




# JAW OPERATION



## Matoset

<b>MB-991-ZEST-590</b>	BP	Tyvek®	pieces
------------------------	----	--------	--------

A -	Drape 200cm×140cm with opening 13.5cm×9.5cm×3.5cm	1
B -	Gauze swab 17-thread 12-ply 5cm×5cm	2
C -	Gauze swab 17-thread 12-ply 10cm×10cm	2
D -	Plastic container 500ml	1
E -	Drain with Yankauer suction tip CH10 300cm	1
F -	Syringe with curved tip 12ml	1
G -	Scalpel blade no. 15	1




## Matoset

<b>MB-991-ZOL3-003*</b>	SMS	PFP	pieces
-------------------------	-----	-----	--------

A -	Drape 160cm×90cm with opening 10cm×15cm and Velcro tape**	1
B -	Drape 90cm×90cm (for set wrapping)	1
C -	Lamp handle cover	1
D -	Gauze swab 17-thread 12-ply 5cm×5cm	5
E -	Drain with Yankauer suction tip 185cm	1
F -	Plastic container 125ml	1
G -	Syringe with curved tip 12ml	1
H -	Scalpel blade no. 15	1

\* information on the laminate given in the code is not consistent with the instructions on page 2; in this case, the main drape is made of SMS polypropylene nonwoven  
 \*\* the Velcro allows you to detach it multiple times if necessary




# HEAD/NECK SURGERY



## Matodrape

MB-135-ZOL2-117	Blue Special	PFP	pieces
-----------------	--------------	-----	--------

A -	Drape 240cm×180cm with opening 15cm×10cm and adhesive tape	1
B -	Instrument table cover 190cm×150cm	1
C -	Mayo table cover 145cm×80cm	1
D -	Medical fixing tape 9cm×50cm	1





# STRUMECTOMY

## Matodrape

MB-135-ZOL2-020	Blue Special	PFP	pieces
-----------------	--------------	-----	--------


A -	Drape 240cm×180cm with adhesive tape	1
B -	Drape 240cm×150cm with U cut-out 63cm×7cm and adhesive tape	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Medical fixing tape 9cm×50cm	1
F -	Nonwoven hand towel 40cm×20cm	2



## Matodrape

MB-135-ZOL3-020	Blue Comfort	PFP	pieces
-----------------	--------------	-----	--------

A -	Drape 240cm×180cm with adhesive tape	1
B -	Drape 240cm×150cm with U cut-out 63cm×7cm and adhesive tape	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Medical fixing tape 9cm×50cm	1
F -	Nonwoven hand towel 40cm×20cm	4





# LARYNGOLOGY

## Matodrape

<b>MB-135-ZOTF-010</b>	TF	PPF	pieces
------------------------	----	-----	--------


- A - Drape 240cm×180cm with opening 10cm and adhesive tape 1
- B - Drape 180cm×150cm (for set wrapping) 1
- C - Nonwoven hand towel 40cm×20cm 1



## Matodrape

<b>MB-135-ZOL2-120</b>	Blue Special	PPF	pieces
------------------------	--------------	-----	--------



- A - Drape 210cm×150cm with opening 10cm 1
- B - Mayo table cover 145cm×80cm 1
- C - Medical fixing tape 9cm×50cm 1



## Matodrape

<b>MB-135-ZOL2-085</b>	Blue Special	PPF	pieces
------------------------	--------------	-----	--------

- A - Drape 210cm×150cm with opening 8cm×6cm and adhesive tape 1
- B - Instrument table cover 190cm×150cm 1
- C - Mayo table cover 145cm×80cm 1
- D - Medical fixing tape 9cm×50cm 2



# PACEMAKER IMPLANTATION




plus



## Matoset

	MA-991-ZOL2-654	Blue Special	Tyvek®	pieces
A -	Drape 320cm×240cm with 2 openings 12cm, incise film and absorbent layer			1
B -	Instrument table cover 190cm×150cm			1
C -	Mayo table cover 145cm×80cm			1
D -	Film cover for medical equipment - round with elastic band (22+60cm)×30cm			1
E -	Film cover for medical equipment - rectangular with elastic band (17+71cm)×91cm			1
F -	Medical fixing tape 9cm×50cm			1
G -	Nonwoven hand towel 40cm×20cm			2
H -	Surgical gown Perfect L			1
I -	Surgical gown Perfect XL			1
J -	Gauze swab 17-thread 8-ply 10cm×10cm			50
K -	Plastic container 125ml			1
L -	Sponge stick for operating field disinfection			2
M -	Retractable surgical scalpel no. 10			1




# ANGIOGRAPHY



## Matodrape

	MB-136-ZOL2-009	Blue Special	PFP	pieces
A -	Drape 320cm×180cm with 2 openings 8cm and incise film			1
B -	Drape 90cm×75cm with adhesive tape			2
C -	Drape 90cm×75cm			2
D -	Drape 150cm×90cm (for set wrapping)			1
E -	Nonwoven hand towel 40cm×20cm			2
F -	Plastic sponge forceps 24cm*			1
G -	Plastic container 250ml			3

\* perfect for holding dressing material




# ANGIOGRAPHY

plus



## Matoset

	MA-991-ZOL2-701	Blue Special	Tyvek®	pieces
A -	Drape 340cm×266cm with 2 openings 12cm, 2 openings 10cm, incise film in all openings, absorbent layer and 2 film panels			1
B -	Drape 75cm×45cm with adhesive tape			1
C -	Drape 200cm×150cm (for set wrapping)			1
D -	Film cover for medical equipment - rectangular with elastic band (17÷71cm)×91cm			1
E -	Film cover for medical equipment - rectangular with elastic band 127×91cm			1
F -	Film cover for medical device cables 15cm×250cm			1
G -	Nonwoven hand towel 40cm×20cm			2
H -	Surgical gown Perfect L			1
I -	Surgical gown Perfect XL			1
J -	Gauze swab 17-thread 8-ply 10cm×10cm			100
K -	Elastic bandage 15cm×5m			5
L -	Needle 1.2mm×40mm			1
M -	Needle 0.7mm×30mm			1
N -	Metal Pean forceps 14cm			1
O -	Plastic container 2500ml for cables			1
P -	Plastic container 500ml			1
Q -	Plastic container 500ml			1
R -	Syringe 2-piece L 20ml			1
S -	Syringe 3-piece L 20ml			1
T -	Syringe 2-piece L 10ml			1
U -	Syringe 3-piece L 10ml			1
V -	Syringe 3-piece LL 10ml			1
W -	Syringe 3-piece LL 5ml			1
X -	Surgical scalpel no. 11			1
Y -	Ultrasound gel in a sachet 20ml			1




plus



Matodrape

	MA-135-ZOL2-438	Blue Special	PFP	pieces
A -	Drape 340cm×240cm with 2 openings 12cm, adhesive tapes, absorbent layer and film panel			1
B -	Drape 90cm×75cm			4
C -	Instrument table cover 190cm×150cm			1
D -	Film cover for medical equipment - rectangular with elastic band (17÷71cm)×91cm			2
E -	Film cover for medical equipment - rectangular with elastic band (28÷102cm)×91cm			1
F -	1-chamber film pouch 40cm×30cm with stiffening element and adhesive tape			1
G -	Surgical gown Comfort Plus XL			2
H -	Nonwoven swab 4-ply 40g/m <sup>2</sup> 10cm×20cm			20
I -	Plastic container 200ml			2
J -	Plastic container 125ml			1
K -	Sponge stick for operating field disinfection			1




plus



Matodrape

	MB-136-ZOL2-016	Blue Special	PFP	pieces
A -	Drape 320cm×235cm with 2 openings 12cm, absorbent layer, incise film and film panel*			1
B -	Drape 150cm×90cm			1
C -	Drape 200cm×180cm (for set wrapping)			1
D -	Film cover for medical equipment - rectangular with elastic band (28÷102cm)×91cm			1
E -	Film cover for medical equipment - round with elastic band (33÷112cm)×56cm			1
F -	Soft Basic underpad 60cm×60cm			1
G -	Surgical gown Standard L			2

\* the transparent film panel protects but does not obstruct the field of vision




# ANGIOGRAPHY

## Matodrape

**MA-135-ZOL3-046** Blue Comfort PFP pieces

A -	Drape 320cm×180cm with 2 openings 8cm and adhesive tapes	1
B -	Drape 200cm×180cm (for set wrapping)	1
C -	Film cover for medical equipment - round with elastic band (33±112cm)×56cm	2
D -	Soft Basic underpad 60cm×60cm	1
E -	Nonwoven hand towel 40cm×20cm	2
F -	Surgical gown Standard L	1
G -	Surgical gown Perfect XL	1
H -	Plastic container 500ml	1
I -	Sponge stick for operating field disinfection	2



plus


## Matoset

**MB-991-ZOL2-042** Blue Special PFP pieces

A -	Drape 330cm×235cm with 2 openings 12cm, adhesive tapes, absorbent layer and film panel*	1
B -	Drape 150cm×150cm (for set wrapping)	1
C -	Film cover for medical equipment - round with elastic band (33±112cm)×56cm	1
D -	Film cover for medical equipment - rectangular with elastic band (28±102cm)×91cm	1
E -	Soft Basic underpad 90cm×60cm	1
F -	Cellulose towel 50cm×40cm	1
G -	Surgical gown Standard L	1
H -	Surgical gown Perfect L	1
I -	Gauze swab 13-thread 8-ply 7.5cm×7.5cm	100
J -	Needle 0.7mm×30mm	1
K -	Metal Pean forceps 14cm	1
L -	Plastic container 2500ml for cables**	1
M -	Plastic container 1000ml	1
N -	Plastic container 500ml	1
O -	Syringe 3-piece LL 20ml	3
P -	Syringe 3-piece LL 10ml	5
Q -	Syringe 3-piece LL 5ml	1
R -	Syringe 3-piece LL 3ml	1
S -	Sponge stick for operating field disinfection	1
T -	Surgical scalpel no. 11	1



metal Pean forceps




\* the transparent film panel protects but does not obstruct the field of vision  
 \*\* this type of container enables cables to be sorted during the procedure


# VARIX SURGERY



## Matodrape



MB-135-ZOL2-045	Blue Special	PFP	pieces
A -	Drape 240cm×180cm with U cut-out 63cm×7cm and adhesive tape		1
B -	Drape 240cm×180cm with adhesive tape		1
C -	Drape 200cm×180cm		1
D -	Instrument table cover 190cm×150cm		1
E -	Mayo table cover 145cm×80cm		1
F -	Medical fixing tape 9cm×50cm		1
G -	Nonwoven hand towel 40cm×20cm		2




## Matoset




MB-991-ZOL2-014	Blue Special	PFP	pieces
A -	Drape 240cm×180cm with U cut-out 63cm×7cm and adhesive tape		1
B -	Drape 240cm×180cm with adhesive tape		1
C -	Drape 200cm×180cm		1
D -	Instrument table cover 190cm×150cm		1
E -	Mayo table cover 145cm×80cm		1
F -	Medical fixing tape 9cm×50cm		1
G -	Nonwoven hand towel 40cm×20cm		2
H -	Gauze swab 17-thread 8-ply 10cm×10cm*		100
I -	Gauze tampon - bean 17-thread with X-ray detectable thread 9.5cm×9.5cm**		5
J -	Gauze abdominal swab 17-thread 6-ply with X-ray detectable thread and cotton loop 45cm×45cm		6

\* tying in stacks by 25 pcs facilitates counting before the procedure




## VARIIX SURGERY



### Matoset

**MB-991-ZOL3-007** Blue Comfort Tyvek® pieces

A -	Drape 240cm×200cm with U cut-out 60cm×7cm and adhesive tape	1
B -	Drape 240cm×180cm with adhesive tape	1
C -	Instrument table cover 190cm×150cm	1
D -	Sheet of swab nonwoven 80cm×60cm	5
E -	Soft Super underpad 90cm×60cm	1
F -	Nonwoven hand towel 80cm×60cm	2
G -	Surgical gown Standard L	2
H -	Nonwoven foot cover with ties	1
I -	Gauze swab 17-thread 16-ply with X-ray detectable thread 7.5cm×7.5cm	50
J -	Elastic bandage 12cm×5m with a clip	1
K -	Elastic bandage 15cm×5m with a clip	3



## CARDIAC SURGERY


plus

### Matodrape

**MB-136-ZOL3-016** Blue Comfort PFP pieces

A -	Drape 325cm×250cm with 2 openings 14cm, incise film, absorbent layer, 1 pouch and film panel*	1
B -	Drape 150cm×120cm (for set wrapping)	1

\* the transparent film panel protects but does not obstruct the field of vision







## UNIVERSAL SETS



## Matodrape

MB-135-ZOTF-007	TF	PFP	pieces
-----------------	----	-----	--------





A -	Drape 300cmx180cm with adhesive tape		3
B -	Drape 90cmx75cm with adhesive tape		2
C -	Mayo table cover 145cmx80cm		1
D -	Medical fixing tape 9cmx50cm		1
E -	Nonwoven hand towel 40cmx20cm		2



## Matodrape

MB-135-ZOTF-011	TF	BLI	pieces
-----------------	----	-----	--------

A -	Drape 200cmx180cm with adhesive tape		1
B -	Drape 180cmx180cm with adhesive tape		1
C -	Drape 90cmx75cm with adhesive tape		2
D -	Instrument table cover 190cmx150cm		1
E -	Mayo table cover 145cmx80cm		1
F -	Medical fixing tape 9cmx50cm		1




# UNIVERSAL SETS

## Matodrape



MB-135-ZOL2-024	Blue Special	BLI	pieces
-----------------	--------------	-----	--------

A -	Drape 240cm×150cm with adhesive tape		1
B -	Drape 180cm×170cm with adhesive tape		1
C -	Drape 90cm×75cm with adhesive tape		2
D -	Nonwoven hand towel 40cm×20cm		2




## Matodrape




MB-135-ZOL2-025	Blue Special	BLI	pieces
-----------------	--------------	-----	--------

A -	Drape 240cm×150cm with adhesive tape		1
B -	Drape 180cm×170cm with adhesive tape		1
C -	Drape 90cm×75cm with adhesive tape		2
D -	Instrument table cover 190cm×150cm		1
E -	Mayo table cover 145cm×80cm		1
F -	Medical fixing tape 9cm×50cm		1
G -	Nonwoven hand towel 40cm×20cm		2




plus

## Matodrape



MB-135-ZOL2-353	Blue Special	BLI	pieces
-----------------	--------------	-----	--------

A -	Drape 240cm×150cm with adhesive tape and absorbent layer		1
B -	Drape 200cm×170cm with adhesive tape and absorbent layer		1
C -	Drape 90cm×75cm with adhesive tape and absorbent layer		2
D -	Instrument table cover 190cm×150cm		1
E -	Mayo table cover 145cm×80cm		1
F -	Medical fixing tape 9cm×50cm		1
G -	Nonwoven hand towel 40cm×20cm		2



# UNIVERSAL SETS




## Matodrape

**MB-135-ZOL2-082** Blue Special PFP pieces

A -	Drape 240cm×180cm with adhesive tape	1
B -	Drape 180cm×170cm with adhesive tape	1
C -	Drape 90cm×75cm with adhesive tape	2
D -	Instrument table cover 190cm×150cm	1
E -	Mayo table cover 145cm×80cm	1
F -	2-chamber film pouch 42cm×35cm with stiffening element* and adhesive tape	1
G -	Medical fixing tape 9cm×50cm	1
H -	Nonwoven hand towel 40cm×20cm	2
I -	Surgical gown Perfect L	1
J -	Surgical gown Perfect XL	2

\* the stiffening element enables adjustment of the opening




## Matodrape

**MB-135-ZOL2-110** Blue Special BLI pieces

A -	Drape 240cm×150cm with adhesive tape	1
B -	Drape 200cm×180cm with adhesive tape	1
C -	Drape 100cm×75cm with adhesive tape	2
D -	Drape 45cm×45cm with adhesive tape	1
E -	Drape 45cm×45cm	4
F -	Drape 200cm×150cm (for set wrapping)	1
G -	Mayo table cover 145cm×80cm	1
H -	1-chamber film pouch 40cm×30cm with stiffening element and adhesive tape	1
I -	Medical fixing tape 9cm×50cm	1
J -	Self-adhesive tape - Velcro 2cm×25cm*	2
K -	Nonwoven hand towel 30cm×30cm	2
L -	Cable organiser 15cm×10cm with adhesive tape**	2

\* the Velcro allows you to detach it multiple times if necessary  
\*\* the organiser helps to organise and hold cables and drains




# UNIVERSAL SETS

## Matodrape

**MB-135-ZOL3-028** Blue Comfort BLI pieces

A -	Drape 240cm×150cm with adhesive tape	1
B -	Drape 180cm×170cm with adhesive tape	1
C -	Drape 90cm×75cm with adhesive tape	2
D -	Instrument table cover 190cm×150cm	1
E -	Mayo table cover 145cm×80cm	1
F -	Medical fixing tape 9cm×50cm	1
G -	Nonwoven hand towel 40cm×20cm	4




## Matoset

**MB-991-ZOL2-072** Blue Special PFP pieces

A -	Drape 240cm×180cm with adhesive tape	1
B -	Drape 180cm×150cm with adhesive tape	1
C -	Drape 180cm×150cm	1
D -	Drape 90cm×75cm with adhesive tape	2
E -	Instrument table cover 190cm×150cm	1
F -	Mayo table cover 145cm×80cm	1
G -	Nonwoven hand towel 40cm×20cm	4
H -	Surgical gown Standard L with 2 hand towels 40cm×20cm*	1
I -	Surgical gown Standard XL	2
J -	Gauze swab 17-thread 8-ply with X-ray detectable thread 7.5cm×7.5cm	20
K -	Gauze tampon - ball 17-thread with X-ray detectable thread 15cm×15cm	40
L -	Syringe 3-piece LL 20ml	1

\* 2 hand towels packed separately allow to separate the clothing from the materials for performing the procedure




# LAPAROTOMY



## Matodrape

<b>MB-135-ZOL2-009</b>	Blue Special	PPF	pieces
------------------------	--------------	-----	--------


- A - Drape 300cm×240cm with opening 10cm×20cm and adhesive tape 1
- B - Nonwoven hand towel 40cm×20cm 2



## Matodrape

<b>MB-135-ZOL2-010</b>	Blue Special	BLI	pieces
------------------------	--------------	-----	--------

- A - Drape 300cm×240cm with opening 10cm×20cm and adhesive tape 1
- B - Instrument table cover 190cm×150cm 1
- C - Mayo table cover 145cm×80cm 1
- D - Medical fixing tape 9cm×50cm 1
- E - Nonwoven hand towel 40cm×20cm 2




# LAPAROSCOPY

## Matodrape

<b>MB-135-ZOL2-046</b>	Blue Special	BLI	pieces
------------------------	--------------	-----	--------


- A - Drape 300cm×180cm with opening 28cm×19cm×20cm and adhesive tape 1
- B - Instrument table cover 190cm×150cm 1
- C - Mayo table cover 145cm×80cm 1
- D - 2-chamber film pouch 42cm×35cm with stiffening element\* and adhesive tape 1
- E - Medical fixing tape 9cm×50cm 1
- F - Nonwoven hand towel 40cm×20cm 2

\* stiffening element enables adjustment of the opening




# LAPAROSCOPY

plus



## Matodrape


	MA-135-ZOL2-437	Blue Special	PFP	pieces
A -	Drape 300cm×230cm with opening 34cm×30cm, adhesive tape and absorbent layer			1
B -	Instrument table cover 190cm×150cm			1
C -	Mayo table cover 145cm×80cm			1
D -	Medical fixing tape 9cm×50cm			1
E -	Cable organiser 15cm×10cm with adhesive tape			2
F -	2-chamber film pouch 42cm×35cm with stiffening element and adhesive tape			1
G -	Nonwoven hand towel 40cm×20cm			2



## Matodrape


	MB-135-ZOL2-048	Blue Special	PFP	pieces
A -	Drape 260cm×240cm with opening 34cm×30cm, adhesive tape and integrated legs*			1
B -	Instrument table cover 190cm×150cm			1
C -	Mayo table cover 145cm×80cm			1
D -	2-chamber film pouch 42cm×35cm with stiffening element and adhesive tape			2
E -	Medical fixing tape 9cm×50cm			1
F -	Nonwoven hand towel 40cm×20cm			2

\* integrated legs = one complete drape instead of 1 main drape and 2 leg covers



## Matodrape

	MB-135-ZOL3-038	Blue Comfort	PFP	pieces
A -	Drape 300cm×180cm with opening 28cm×19cm×20cm and adhesive tape			1
B -	Instrument table cover 190cm×150cm			1
C -	Mayo table cover 145cm×80cm			1
D -	1-chamber film pouch 40cm×30cm with stiffening element and adhesive tape			1
E -	Medical fixing tape 9cm×50cm			1
F -	Nonwoven hand towel 40cm×20cm			4




# THORACOTOMY



## Matodrape

**MB-135-ZOL2-042** Blue Special PFP pieces


A -	Drape 300cm×180cm with cut-out 60cm×30cm and adhesive tape	1
B -	Drape 180cm×90cm with adhesive tape	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Nonwoven hand towel 40cm×20cm	2



## Matodrape

**MB-135-ZOL3-047** Blue Comfort PFP pieces

A -	Drape 240cm×200cm with adhesive tape	1
B -	Drape 240cm×180cm with adhesive tape	1
C -	Drape 90cm×75cm with adhesive tape	2
D -	Drape 75cm×45cm	4
E -	Drape SMS 90cm×80cm	2
F -	Instrument table cover 190cm×150cm	1
G -	Mayo table cover 145cm×80cm	2
H -	2-chamber film pouch 42cm×35cm with stiffening element and adhesive tape	1
I -	Surgical gown Perfect M	1
J -	Surgical gown Perfect L	1




## Matodrape



MB-135-ZOL2-012	Blue Special	PFP	pieces
-----------------	--------------	-----	--------

A -	Drape 200cm×180cm with opening 20cm×10cm and adhesive tape	1
B -	Drape 90cm×75cm	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Limb cover 120cm×75cm	2
F -	Nonwoven hand towel 40cm×20cm	2




## Matodrape



MB-135-ZOL3-011	Blue Comfort	PFP	pieces
-----------------	--------------	-----	--------

A -	Drape 200cm×180cm with opening 20cm×10cm and adhesive tape	1
B -	Drape 90cm×75cm	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Limb cover 120cm×75cm	2
F -	Nonwoven hand towel 40cm×20cm	4




## Matodrape set for laparoscopic gynaecological operations



MB-135-ZOL2-040	Blue Special	PFP	pieces
-----------------	--------------	-----	--------

A -	Drape 265cm×180cm with 2 openings 28cm×20cm×19cm and 15cm×10cm (the bottom opening is additionally covered with a flap) and <b>integrated legs*</b>	1
B -	Drape 90cm×75cm	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Medical fixing tape 9cm×50cm	1
F -	Nonwoven hand towel 40cm×20cm	2

\* integrated legs = one complete drape instead of 1 main drape and 2 leg covers




**Matodrape**

**MB-135-ZOL3-034** Blue Comfort BLI pieces



A -	Drape 245cm×180cm with opening 15cm×10cm, adhesive tape and integrated legs*	1
B -	Drape 90cm×75cm	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Medical fixing tape 9cm×50cm	2
F -	Nonwoven hand towel 40cm×20cm	4

\* integrated legs = one complete drape instead of 1 main drape and 2 leg covers




**Matoset**

**MB-991-ZOL2-172** Blue Special PFP pieces



A -	Drape 275cm×200cm with opening 10cm×15cm, adhesive tape, integrated legs and liquid collection pouch with a sieve	1
B -	Instrument table cover 190cm×150cm	1
C -	Mayo table cover 145cm×80cm	1
D -	Medical fixing tape 9cm×50cm	1
E -	Soft Super underpad 90cm×60cm	1
F -	Surgical gown Perfect L	1
G -	Surgical gown Perfect XL	2
H -	Gauze swab 17-thread 12-ply with X-ray detectable thread 10cm×20cm	5
I -	Gauze swab 17-thread 12-ply with X-ray detectable thread 7.5cm×7.5cm	50
J -	Gauze tampon - ball 17-thread with X-ray detectable thread 15cm×15cm	20
K -	Gauze abdominal swab 17-thread 4-ply with X-ray detectable thread and cotton loop 45cm×45cm	5




# GYNAECOLOGY



## Matoset

**MB-991-ZOL3-103** Blue Comfort Tyvek® pieces

A -	Drape 250cm×250cm with opening 12cm×8cm, adhesive tape integrated legs	1
B -	Instrument table cover 190cm×150cm	1
C -	Mayo table cover 145cm×80cm	2
D -	Incise film 10cm×20cm	1
E -	Gauze swab 17-thread 12-ply with X-ray detectable thread 10cm×10cm	10
F -	Gauze swab 17-thread 12-ply with X-ray detectable thread 7.5cm×7.5cm	40
G -	Gauze ribbon 17-thread 4ply with X-ray detectable thread 5cm×2m	2
H -	Monopolar cautery	1
I -	Rigid blister container 10cm×20cm×4cm	1
J -	Syringe 3-part L	1
K -	Scalpel blade no. 20	2




## Matoset

**MB-991-ZOL2-164** Blue Special Tyvek® pieces

A -	Drape 275cm×200cm with opening 10cm×15cm, adhesive tape, integrated legs and liquid collection pouch with a sieve	1
B -	Drape 90cm×75cm	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Film cover for medical device cables 15cm×250cm	1
F -	Nonwoven hand towel 40cm×20cm	2
G -	Surgical gown Perfect L	1
H -	Surgical gown Perfect XL	1
I -	Gauze swab 17-thread 12-ply with X-ray detectable thread 10cm×10cm*	20
J -	Gauze tampon - ball 17-thread 50cm×50cm	3
K -	Rigid blister container 10cm×20cm×4cm	1

\* tying in stacks by 10 pcs facilitates counting before the procedure



# ABDOMINAL-PERINEAL SURGERY




## Matodrape

**MB-136-ZOL3-008** Blue Comfort PFP pieces



A -	Drape 260cm×240cm with 2 openings (upper opening 34cm×30cm filled with incise film, the bottom opening 15cm×10cm covered with a flap) and integrated leg*	1
B -	Drape 90cm×75cm	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	2-chamber film pouch 42cm×35cm with stiffening element and adhesive tape	1
F -	Medical fixing tape 9cm×50cm	1
G -	Nonwoven hand towel 40cm×20cm	4

\* integrated legs = one complete drape instead of 1 main drape and 2 leg covers




## DELIVERY

## Matodrape

**MB-135-ZOMI-011** MIX PFP pieces

A -	Drape L2 130cm×90cm	1
B -	Drape L2 150cm×90cm (for set wrapping)	1
C -	Nonwoven towel 80cm×60cm	2
D -	Soft Super underpad 90cm×60cm	1
E -	Soft Super underpad 60cm×60cm	1
F -	Flannel blanket 160cm×75cm*	1
G -	Baby hat	1

\* a variety of designs




# DELIVERY

## Matodrape

MB-135-ZOMI-008	MIX	PFP	pieces
-----------------	-----	-----	--------


A -	Under buttocks drape L2 113cm×90cm liquid collection pouch of laminate and stiffening element	1
B -	Drape L3 150cm×90cm (for set wrapping)	1
C -	Nonwoven towel 80cm×60cm	1
D -	Nonwoven hand towel 25cm×20cm	6
E -	Soft Super underpad 90cm×60cm	1
F -	Soft Super underpad 60cm×60cm	1



## Matoset

MB-991-ZBS2-002	Blue Special	Tyvek®	pieces
-----------------	--------------	--------	--------

A -	Drape L2 150cm×90cm (for set wrapping)	1
B -	Limb cover 120cm×80cm	1
C -	Nonwoven towel 80cm×60cm	1
D -	Soft Super underpad 90cm×60cm	1
E -	Flannel blanket 160cm×75cm*	1
F -	Baby hat	1
G -	Examination glove powder-free M	2
H -	Nonwoven swab 4-ply 40g/m <sup>2</sup> 7.5cm×7.5cm	10
I -	Nonwoven swab 4-ply 40g/m <sup>2</sup> 10cm×20cm	5
J -	Gauze tampon - ball 17-thread 30cm×30cm	5
K -	Baby measuring tape	1
L -	Umbilical cord clamp	2
M -	Metal episiotomy scissors 18cm**	1
N -	Metal umbilical scissors 10.5cm	1
O -	Plastic sponge forceps 24cm***	1
P -	Plastic kidney dish 700ml	1



\* a variety of designs  
 \*\* metal scissors with a safety ball tip  
 \*\*\* ideal for holding dressing material

## DELIVERY




### Matoset

<b>MB-991-ZOL2-048</b>	Blue Special	PFP	pieces
------------------------	--------------	-----	--------



A - Under buttocks drape L2 113cm×90cm liquid collection pouch of laminate and stiffening element	1
B - Drape TF 90cm×75cm	1
C - Nonwoven hand towel 40cm×35cm	3
D - Soft Super underpad 60cm×60cm	1
E - Surgical gown Perfect L	1
F - Gauze swab 17-thread 8-ply 7.5cm×7.5cm	10



## CAESAREAN SECTION


### Matodrape

<b>MB-136-ZOL2-003</b>	Blue Special	PFP	pieces
------------------------	--------------	-----	--------



A - Drape 325cm×250cm with opening 35cm×32cm, liquid collection pouch*, incise film and cable organisers**	1
B - Instrument table cover 190cm×150cm	1
C - Mayo table cover 145cm×80cm	1
D - Soft Basic underpad 90cm×60cm	1
E - Nonwoven hand towel 40cm×20cm	2

\* the pouch with valves for draining fluids  
 \*\* the organiser helps to organise and hold cables and drains




# CAESAREAN SECTION

## Matoset



MA-991-ZOL2-449	Blue Special	PFP	pieces
A -	Drape 270cm×200cm with opening 35cm×32cm, liquid collection pouch, incise film and cut-out in the film and with integrated legs		1
B -	Drape 100cm×105cm		1
C -	Instrument table cover 190cm×150cm		1
D -	Mayo table cover 145cm×80cm		1
E -	Medical fixing tape 9cm×50cm		1
F -	Cellulose towel 50cm×40cm		3
G -	Flannel blanket 80cm×75cm		1
H -	Baby hat		1
I -	Surgical gown Standard M		1
J -	Surgical gown Perfect L		1
K -	Surgical gown Perfect XL		1
L -	Gauze swab 17-thread 8-ply with X-ray detectable thread 10cm×10cm*		10
M -	Gauze swab 17-thread 8-ply with X-ray detectable thread 10cm×20cm*		20
N -	Gauze abdominal swab 17-thread 4-ply with X-ray detectable thread and cotton loop 45cm×45cm		4

\* tying in stacks by 10 pcs facilitates counting before the procedure




## Matodrape



MB-136-ZOL2-010	Blue Special	PFP	pieces
A -	Drape 320cm×180cm with opening 35cm×32cm, liquid collection pouch*, incise film and cable organisers**		1
B -	Instrument table cover 190cm×150cm		1
C -	Mayo table cover 145cm×80cm		1
D -	Soft Basic underpad 90cm×60cm		1
E -	Nonwoven hand towel 40cm×20cm		2

\* the pouch with valves for draining fluids  
\*\* the organiser helps to organise and hold cables and drains



# CAESAREAN SECTION




## Matodrape

**MB-136-ZOL2-015** Blue Special PFP pieces



A -	Drape 320cm×180cm with opening 35cm×32cm, liquid collection pouch*, incise film and cable organisers**	1
B -	Instrument table cover 190cm×150cm	1
C -	Mayo table cover 145cm×80cm	1
D -	Soft Super underpad 90cm×60cm	1
E -	Surgical gown Perfect L	1
F -	Surgical gown Perfect XL	2




\* the pouch with valves for draining fluids  
\*\* the organiser helps to organise and hold cables and drains

## Matodrape

**MB-136-ZOL3-002** Blue Comfort PFP pieces



A -	Drape 325cm×250cm with opening 35cm×32cm, liquid collection pouch*, incise film and cable organisers**	1
B -	Instrument table cover 190cm×150cm	1
C -	Mayo table cover 145cm×80cm	1
D -	Medical fixing tape 9cm×50cm	1
E -	Nonwoven hand towel 40cm×20cm	4
F -	Flannel blanket 160cm×75cm***	1
G -	Baby hat	1



\* the pouch with valves for draining fluids  
\*\* the organiser helps to organise and hold cables and drains  
\*\*\* a variety of designs

# BABY SETS

## Matodrape

**MB-135-ZOMI-010** MIX PFP pieces



A -	Soft Super underpad 60cm×60cm	1
B -	Nonwoven towel 80cm×60cm	2
C -	Flannel blanket 160cm×75cm*	1
D -	Baby hat	1





\* a variety of designs

# OPHTHALMOLOGY

## Matoset

MA-991-ZES2-846	Blue Special	PFP	pieces
-----------------	--------------	-----	--------




A -	Drape 90cm×75cm with opening 8cm×6cm and adhesive tape	1
B -	Drape 150cm×90cm (for set wrapping)	1
C -	Nonwoven swab 4-ply 30g/m <sup>2</sup> 10cm×10cm	5
D -	Gauze tampon - ball 17-thread 20cm×20cm	3
E -	Metal barraquer for adults 4cm, closed loops	1
F -	Plastic scleral marker for marking injection sites with 2 tips: 3.5mm-4.0 mm	1
G -	Ocular stick Pro Ophta 5mm×64mm	2



## Matoset

MA-991-ZOTM-013	SMS	PFP	pieces
-----------------	-----	-----	--------

A -	Drape 90cm×80cm with opening 10cm×7cm, <b>incise film</b> and slit in the film	1
B -	Drape 90cm×75cm (for set wrapping)	1
C -	Gauze swab 17-thread 12-ply 7.5cm×7.5cm	5
D -	Gauze eye dressing pad 17-thread 5cm×6cm	1
E -	Metal barraquer for adults 4cm, closed loops	1
F -	Plastic container 125ml	1
G -	Syringe 2-piece L 10ml	1
H -	Syringe 2-piece L 2ml	1
I -	Plastic Pean forceps 13cm	1
J -	Plastic cotton swab 15cm	1
K -	Plastic scleral marker for marking injection sites with 2 tips: 3.5mm-4.0 mm	1






## Matoset

MA-991-ZOTF-023	TF	PPF	pieces
-----------------	----	-----	--------

A -	Drape 180cm×150cm with opening 10cm×7cm with incise film		1
B -	Drape SMS 130cm×80cm		1
C -	Drape 75cm×45cm		1
D -	Drape 130cm×90cm (for set wrapping)		1
E -	Mayo table cover 145cm×80cm		1
F -	Armrest cover 60cm×40cm with elastic band		2
G -	Film cover for medical equipment – round with elastic band 72cm×36cm		1
H -	1-chamber film pouch 35cm×18.5cm with stiffening element and adhesive tape		1
I -	Nonwoven hand towel 40cm×20cm		4
J -	Surgical gown Standard L		2
K -	Nonwoven swab 8-ply 30g/m <sup>2</sup> 7.5cm×7.5cm		10
L -	Nonwoven adhesive dressing with absorbent pad 6.5cm×9.5cm		1
M -	Plastic container 120ml		1
N -	Plastic container 60ml		1
O -	Syringe 3-part LL 5ml		1
P -	Syringe 3-part LL 3ml		2
Q -	Eye spear for ophthalmic procedures 7cm – packaging by 5 pcs		1
R -	Plastic Pean forceps 13cm		1



-- A



-- B



-- C



-- D



-- E



-- F




-- G



-- H



-- I



-- J



10x -- K



-- L



-- M



-- N



-- O



-- P

5x



-- Q



-- R


# OPHTHALMOLOGY



## Matoset

	MB-991-ZOL2-016	Blue Special	BLI	pieces
A -	Drape 150cm×150cm with opening 8cm, incise film and 2 liquid collection pouches*			1
B -	Instrument table cover 190cm×150cm			1
C -	Nonwoven hand towel 80cm×60cm			2
D -	Surgical gown Standard L			2
E -	Gauze swab 17-thread 16-ply 7.5cm×7.5cm			10
F -	Needle 1.2mm×40mm			1
G -	Needle 0.5mm×25mm			1
H -	Plastic eye shield			1
I -	Plastic container 125ml			1
J -	Syringe 3-piece LL 20ml			1
K -	Syringe 3-piece LL 3ml			3
L -	Syringe 3-piece LL 1ml			1


\* the pouches with stiffening element for better opening adjustment



## Matoset

	MA-991-ZOL2-661	Blue Special	PFP	pieces
A -	Drape 150cm×150cm with opening 8cm, incise film and 2 liquid collection pouches*			1
B -	Drape 150cm×120cm (for set wrapping)			1
C -	Armrest cover 60cm×40cm with elastic band			2
D -	Surgical gown Standard M			1
E -	Surgical gown Standard L			1
F -	Gauze swab 17-thread 8-ply 5cm×5cm			12
G -	Gauze eye dressing pad 17-thread 2-ply 5cm×6cm			1
H -	Plastic cup 30ml			1
I -	Syringe 3-piece L 20ml			1
J -	Syringe 3-piece L 5ml			1
K -	Syringe 3-piece LL 3ml			2

\* the pouches with stiffening element for better opening adjustment




# KNEE ARTHROSCOPY



## Matodrape

<b>MB-135-ZOTF-001</b>	TF	PPF	pieces
------------------------	----	-----	--------


A -	Drape 180cm×150cm with semi-circular cut-out and adhesive tapes	1
B -	Drape 180cm×150cm with semi-circular cut-out and adhesive tape	1
C -	Limb cover 80cm×45cm	1
D -	Medical fixing tape 9cm×50cm	2



## Matodrape

<b>MB-135-ZOL2-003</b>	Blue Special	PPF	pieces
------------------------	--------------	-----	--------

A -	Drape 320cm×240cm with opening 6cm	1
B -	Limb cover 75cm×37cm	1
C -	Medical fixing tape 9cm×50cm	2
D -	Nonwoven hand towel 40cm×20cm	2




# KNEE ARTHROSCOPY

plus


## Matodrape

	MA-135-ZOL2-436	Blue Special	PFP	pieces
A -	Drape 320cm×230cm with opening 6cm, elastic element and absorbent layer			1
B -	Instrument table cover 190cm×150cm			1
C -	Mayo table cover 145cm×80cm			1
D -	Limb cover 75cm×37cm			1
E -	Medical fixing tape 9cm×50cm			2
F -	Nonwoven hand towel 40cm×20cm			2



## Matodrape

	MB-135-ZOL2-004	Blue Special	PFP	pieces
A -	Drape 320cm×240cm with opening 6cm, elastic element and liquid collection pouch*			1
B -	Instrument table cover 190cm×150cm			1
C -	Mayo table cover 145cm×80cm			1
D -	Limb cover 75cm×37cm			1
E -	Medical fixing tape 9cm×50cm			2
F -	Nonwoven hand towel 40cm×20cm			2



\* the pouch with a valve for draining fluids

# KNEE ARTHROSCOPY




## Matodrape

<b>MB-135-ZOL3-004</b>	Blue Comfort	PPF	pieces
------------------------	--------------	-----	--------



A -	Drape 320cm×240cm with opening 6cm, elastic element and liquid collection pouch*	1
B -	Instrument table cover 190cm×150cm	1
C -	Mayo table cover 145cm×80cm	1
D -	Limb cover 75cm×37cm	1
E -	Medical fixing tape 9cm×50cm	2
F -	Nonwoven hand towel 40cm×20cm	4

\* the pouch with a valve for draining fluids




## Matoset

<b>MB-991-ZOL2-155</b>	Blue Special	PPF	pieces
------------------------	--------------	-----	--------



A -	Drape 300cm×180cm with U cut-out 100cm×7cm and adhesive tape	1
B -	Drape 240cm×150cm with adhesive tape	1
C -	Instrument table cover 190cm×150cm	1
D -	Limb cover 75cm×45cm	1
E -	Film cover for medical device cables 15cm×250cm*	1
F -	Medical fixing tape 9cm×50cm	2
G -	Cellulose hand towel 50cm×40cm	2
H -	Surgical gown Standard L	2
I -	Gauze swab 13-thread 8-ply 10cm×10cm	20
J -	Elastic bandage 15cm×5m with a clip	1




# KNEE ARTHROSCOPY

## Matoset

MB-991-ZOL2-027	Blue Special	Tyvek®	pieces
A -	Drape 320cm×240cm with opening 6cm and elastic element		1
B -	Instrument table cover 190cm×150cm		1
C -	Limb cover 75cm×45cm		1
D -	Film cover for medical device cables 15cm×250cm		2
E -	Medical fixing tape 9cm×50cm		2
F -	Nonwoven hand towel 40cm×20cm		3
G -	Surgical gown Standard M		1
H -	Surgical gown Perfect XL		2
I -	Gauze swab 17-thread 12-ply with X-ray detectable thread 7.5cm×7.5cm		20
J -	Elastic bandage 15cm×4m with a clip		1
K -	Synthetic orthopaedic underpad 10cm×3m		1
L -	Redon drain CH10 50cm*		1
M -	Redon bottle 200ml*		1
N -	Lumbar puncture needle 1.2mm×90mm		1
O -	Metal needle holder 13cm		1
P -	Plastic sponge forceps 19cm**		2
Q -	Plastic Pean forceps 13cm		3
R -	Metal surgical scissors sharp-sharp 11cm		1
S -	Surgical scalpel no. 11		1



\* for draining fluids  
\*\* ideal for holding dressing material




## Matodrape

**MB-135-ZOL3-017** Blue Comfort PFP pieces



A -	Drape 260cm×240cm with U cut-out 100cm×20cm and adhesive tape	1
B -	Drape 240cm×150cm with adhesive tape	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Limb cover 120cm×37cm	1
F -	1-chamber film pouch 40cm×30cm with stiffening element* and adhesive tape	1
G -	Medical fixing tape 9cm×50cm	6
H -	Nonwoven hand towel 40cm×20cm	4

\* the stiffening element for better opening adjustment




## Matodrape

**MA-135-ZOL2-435** Blue Special PFP pieces




A -	Drape 320cm×240cm with 2 openings 10cm×15cm with 2 elastic elements	1
B -	Drape 150cm×80cm	2
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Limb cover 120cm×37cm	2
F -	Medical fixing tape 9cm×50cm	6




# HIP SURGERY

plus



## Matodrape


	MA-135-ZOL2-439	Blue Special	PFP	pieces
A -	Drape 260cm×230cm with U cut-out 100cm×20cm, adhesive tape and absorbent layer			1
B -	Drape 240cm×150cm with adhesive tape and absorbent layer			1
C -	Instrument table cover 190cm×150cm			1
D -	Mayo table cover 145cm×80cm			1
E -	Limb cover 120cm×37cm			1
F -	1-chamber film pouch 40cm×30cm with stiffening element and adhesive tape			1
G -	Medical fixing tape 9cm×50cm			6
H -	Nonwoven hand towel 40cm×20cm			2



## Matoset

	MA-991-ZOL3-189	Blue Comfort	Tyvek®	pieces
A -	Drape 260cm×240cm with U cut-out 100cm×20cm and adhesive tape			1
B -	Drape 240cm×150cm with adhesive tape			1
C -	Drape 240cm×150cm with adhesive tape			1
D -	Drape 180cm×150cm			2
E -	Instrument table cover 190cm×150cm			1
F -	Limb cover 120cm×37cm			1
G -	1-chamber film pouch 40cm×30cm with stiffening element and adhesive tape			1
H -	Medical fixing tape 9cm×50cm			4
I -	Incise film 55cm×45cm			1
J -	Surgical gown Perfect XL			1
K -	Surgical gown Perfect XXL			1
L -	Gauze swab 17-thread 12-ply 10cm×10cm*			80
M -	Gauze abdominal swab 17-thread 4-ply with X-ray detectable thread 45cm×45cm			6
N -	Sponge for cleaning electrodes			1
O -	Yankauer suction tip CH28			1
P -	Drain for suction of the operating field CH30 300cm			1
Q -	Redon drain CH14 70cm			2
R -	Redon bottle 400ml			2
S -	Syringe 3-piece L100ml			1

\* tying in stacks by 10 pcs facilitates counting before the procedure




## Matoset

MB-991-ZOL3-057	Blue Comfort	Tyvek®	pieces
A -	Drape 260cm×240cm with U cut-out 100cm×20cm and adhesive tape		1
B -	Drape 240cm×180cm with adhesive tape		1
C -	Drape 180cm×150cm		1
D -	Drape 90cm×75cm		3
E -	Instrument table cover 190cm×150cm		1
F -	Mayo table cover 145cm×80cm		1
G -	Limb cover 120cm×37cm		1
H -	1-chamber film pouch 40cm×30cm with stiffening element* and adhesive tape		1
I -	Medical fixing tape 9cm×50cm		2
J -	Nonwoven hand towel 40cm×20cm		6
K -	Surgical gown Perfect L with 2 hand towels 40cm×20cm**		1
L -	Surgical gown Perfect XL		2
M -	Surgical gown Perfect XXL		1
N -	Gauze swab 17-thread 8-ply with X-ray detectable thread 10cm×10cm		30
O -	Gauze tampon - ball 17-thread with X-ray detectable thread 30cm×30cm		60
P -	Gauze abdominal swab 17-thread 4-ply with X-ray detectable thread 45cm×45cm		10
Q -	Yankauer suction tip CH28		1
R -	Redon drain CH18 75cm		1
S -	Drain for suction of the operating field CH28 350cm		1
T -	Redon bottle 600ml with connecting drain		1
U -	Syringe 3-piece L 100ml		1
V -	Scalpel blade no. 22		3



\* the stiffening element for better opening adjustment  
 \*\* 2 hand towels packed separately allow to separate the clothing from the materials for performing the procedure





## VERTICAL ISOLATION SETS

### Matodrape

MA-136-ZOPE-001	PE	PFP	pieces
-----------------	----	-----	--------



A - Vertical isolation drape 330cm×240cm with opening 25cm×86cm, incise film, 2 pouches, cable organisers and liquid collection pouch			1
B - Instrument table cover 190cm×150cm			1
C - Mayo table cover 145cm×80cm			1
D - Nonwoven hand towel 40cm×20cm			2



### Matodrape

MA-136-ZOPE-002	PE	PFP	pieces
-----------------	----	-----	--------

A - Vertical isolation drape 330cm×240cm with opening 30cm×40cm, incise film, 2 pouches, cable organisers and liquid collection pouch			1
B - Instrument table cover 190cm×150cm			1
C - Mayo table cover 145cm×80cm			1
D - Nonwoven hand towel 40cm×20cm			2



# SHOULDER SURGERY




## Matodrape



<b>MB-135-ZOL2-049</b>	Blue Special	PPF	pieces
------------------------	--------------	-----	--------

A -	Drape 240cm×200cm with U cut-out 60cm×7cm and adhesive tape	1
B -	Drape 240cm×150cm with adhesive tape	1
C -	Drape 180cm×150cm	1
D -	Instrument table cover 190cm×150cm	1
E -	Mayo table cover 145cm×80cm	1
F -	Limb cover 75cm×37cm	1
G -	Medical fixing tape 9cm×50cm	2
H -	Nonwoven hand towel 40cm×20cm	2




## Matodrape



<b>MB-135-ZOL3-035</b>	Blue Comfort	PPF	pieces
------------------------	--------------	-----	--------

A -	Drape 320cm×240cm with opening 6cm, elastic element and liquid collection pouch*	1
B -	Drape 90cm×75cm with adhesive tape	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Limb cover 75cm×37cm	1
F -	Medical fixing tape 9cm×50cm	2
G -	Nonwoven hand towel 40cm×20cm	4

\* the pouch with a valve for draining fluids




## UPPER LIMB SURGERY

### Matodrape

MB-135-ZOL2-051	Blue Special	PFP	pieces
-----------------	--------------	-----	--------

A -	Drape 260cm×240cm with opening 5cm, elastic element and attached drape 150cm×90cm	1
B -	Instrument table cover 190cm×150cm	1
C -	Mayo table cover 145cm×80cm	1
D -	Limb cover 75cm×37cm	1
E -	Medical fixing tape 9cm×50cm	1
F -	Nonwoven hand towel 40cm×20cm	2




## HAND/FOOT SURGERY

### Matodrape

MB-135-ZOTF-014	TF	PFP	pieces
-----------------	----	-----	--------


A -	Drape 200cm×180cm with opening 4cm and elastic element	1
B -	Drape 180cm×150cm (for set wrapping)	1
C -	Mayo table cover 145cm×80cm	1
D -	Surgical gown Standard L	1
E -	Surgical gown Standard XL	1



### Matodrape

MB-135-ZOL3-030	Blue Comfort	BLI	pieces
-----------------	--------------	-----	--------

A -	Drape 300cm×240cm with opening 4cm and elastic element	1
B -	Instrument table cover 190cm×150cm	1
C -	Mayo table cover 145cm×80cm	1
D -	Nonwoven hand towel 40cm×20cm	4




Matoset

<b>MB-991-ZOL2-026</b>	Blue Special	PFP	pieces
------------------------	--------------	-----	--------



A - Under buttocks drape L2 113cm×90cm liquid collection pouch of laminate and stiffening element*	1
B - Limb cover 130cm×40cm	2
C - Surgical gown Standard L	1
D - Gauze swab 17-thread 12-ply 7.5cm×7.5cm	10

\* the stiffening element for better opening adjustment




Matoset

<b>MB-991-ZOL2-030</b>	Blue Special	PFP	pieces
------------------------	--------------	-----	--------



A - Drape 200cm×180cm with opening 10cm×20cm and adhesive tape	1
B - Drape 150cm×90cm (for set wrapping)	1
C - Limb cover 120cm×80cm	2
D - Gauze swab 17-thread 12-ply 10cm×10cm	5




Matoset

<b>MA-991-ZES2-439</b>	Blue Special	PFP	pieces
------------------------	--------------	-----	--------




A - Drape 275cm×200cm with opening 10cm×15cm, adhesive tape, integrated legs and liquid collection pouch with a sieve	1
B - Instrument table cover 190cm×150cm	1
C - Gauze swab 17-thread 12-ply 10cm×10cm	20
D - Gauze ribbon 17-thread 4-ply 5cm×2m	1
E - Dallop NM urological tape 45cm	1
F - Syringe 2-piece L 20ml	1
G - Scalpel blade no. 20	1




## UROLOGIC SURGERY

### Matoset

MA-991-ZES2-923	Blue Special	PFP	pieces
-----------------	--------------	-----	--------



A -	Drape 200cm×180cm with opening 10cm×20cm and adhesive tape		1
B -	Instrument table cover 190cm×150cm		1
C -	Mayo table cover 145cm×80cm		1
D -	Limb cover 120cm×80cm		2
E -	Medical fixing tape 9cm×50cm		2
F -	Soft Super underpad 90cm×60cm		1
G -	Gauze swab 17-thread 16-ply with X-ray detectable thread 10cm×10cm		10
H -	Gauze tampon - ball 17-thread 30cm×30cm		5
I -	Gauze ribbon 17-thread 4-ply with X-ray detectable thread 2cm×2m		1
J -	Gauze ribbon 17-thread 4-ply with X-ray detectable thread 5cm×1m		1
K -	Dallop NM urological tape 45cm		1
L -	Plastic sponge forceps 24cm		1



## TUR UROLOGIC SURGERY


### Matodrape

MB-135-ZOL2-031	Blue Special	PFP	pieces
-----------------	--------------	-----	--------



A -	Drape 240cm×100cm with 2 openings 5cm and 4cm, semi-circular cut-outs, elastic finger protector and liquid collection pouch with a sieve*		1
B -	Drape 90cm×75cm		1
C -	Instrument table cover 190cm×150cm		1
D -	Mayo table cover 145cm×80cm		1
E -	Limb cover 120cm×75cm		2
F -	Medical fixing tape 9cm×50cm		4
G -	Nonwoven hand towel 40cm×20cm		2

\* the pouch with a sieve for the removed residues and with a valve for draining excess fluids




Matodrape

<b>MB-135-ZOL2-065</b>	Blue Special	PFP	pieces
------------------------	--------------	-----	--------



A -	Drape 215cm×180cm with 2 openings 5cm and 4cm, elastic finger protector, liquid collection pouch with a sieve* and integrated legs**	1
B -	Instrument table cover 190cm×150cm	1
C -	Mayo table cover 145cm×80cm	1
D -	Medical fixing tape 9cm×50cm	1
E -	Nonwoven hand towel 40cm×20cm	2

\* the pouch with a sieve for the removed residues and with a valve for draining excess fluids  
 \*\* integrated legs = one complete drape instead of 1 main drape and 2 leg covers




Matoset

<b>MB-991-ZOTF-003</b>	TF	PFP	pieces
------------------------	----	-----	--------




A -	Drape 180cm×75cm with opening 10cm×20cm, semi-circular cut-outs and adhesive tape	1
B -	Drape 130cm×90cm (for set wrapping)	1
C -	Limb cover 120cm×80cm	2
D -	Surgical gown Standard XL	1
E -	Gauze swab 17-thread 8-ply 5cm×5cm	10
F -	Plastic sponge forceps 24cm*	1

\* ideal for holding dressing material



# PCNL UROLOGIC SURGERY


## Matodrape



MB-135-ZOL3-023	Blue Comfort	PFP	pieces
-----------------	--------------	-----	--------

- |  |   |
|--|---|
| A - Drape 300cm×180cm with opening 22cm×15cm, adhesive tape and liquid collection pouch* | 1 |
| B - Instrument table cover 190cm×150cm   | 1 |

\* the stiffening element for better opening adjustment



# NEPHROSTOMY


## Matoset



MB-991-ZOL2-071	Blue Special	PFP	pieces
-----------------	--------------	-----	--------

- |   |    |
|---|----|
| A - Drape 180cm×170cm with opening 22cm×15cm and adhesive tape                | 1  |
| B - Instrument table cover 190cm×150cm  | 1  |
| C - Film cover for medical device cables 15cm×250cm                           | 1  |
| D - 1-chamber film pouch 40cm×30cm with stiffening element* and adhesive tape | 1  |
| E - Nonwoven hand towel 40cm×20cm   | 2  |
| F - Surgical gown Standard L  | 1  |
| G - Surgical gown Standard XL   | 1  |
| H - Surgical glove powder-free 6.5  | 2  |
| I - Gauze swab 13-thread 12-ply 7.5cm×7.5cm                                   | 20 |

\* the stiffening element for better opening adjustment



# PROSTATE ELECTRORESECTION




## Matodrape



	MB-135-ZOMI-003	TF	PFP	pieces
A -	Drape TF 240cm×180cm with opening 8cm, adhesive tapes* and slits			1
B -	Surgical gown Perfect XL			1
C -	Funnel pouch**			1

\* the adhesive tapes allow the drape to be attached around the patient's legs to form the legs  
 \*\* for attaching to the front of the operator for draining fluids and remains during the procedure




# CYSTOSCOPY

## Matodrape



	MB-135-ZOTF-002	TF	BLI	pieces
A -	Drape 180cm×75cm with opening 20cm×10cm, with semi-circular cut-outs and adhesive tape			1
B -	Instrument table cover 190cm×150cm			1
C -	Limb cover 120cm×60cm			2
D -	Nonwoven hand towel 40cm×20cm			2




# CYSTOSCOPY

## Matodrape

<b>MB-135-ZOL2-008</b>	Blue Special	BLI	pieces
------------------------	--------------	-----	--------



A -	Drape 180cm×75cm with opening 20cm×10cm, semi-circular cut-outs and adhesive tape	1
B -	Drape 90cm×75cm	1
C -	Instrument table cover 190cm×150cm	1
D -	Mayo table cover 145cm×80cm	1
E -	Limb cover 120cm×75cm	2
F -	Nonwoven hand towel 40cm×20cm	2




## Matodrape

<b>MB-135-ZOL3-006</b>	Blue Comfort	PPF	pieces
------------------------	--------------	-----	--------



A -	Drape 180cm×75cm with opening 20cm×10cm, semi-circular cut-outs and adhesive tape	1
B -	Drape 90cm×75cm	1
C -	Limb cover 120cm×75cm	2
D -	Nonwoven hand towel 40cm×20cm	4



Matodress



DISPOSABLE

- compliance with EN 13795 standard
- safe and functional design
- high barrier to liquid penetration
- high absorption capacity in the critical product area
- high strength parameters

Disposable surgical gowns - Matodress

	Standard performance				High performance			
name	Standard	Standard PLUS*	Comfort	Comfort PLUS*	Perfect	Perfect PLUS*	Comfort HP	Comfort HP PLUS*
raw material	polypropylene nonwoven		cotton-like nonwoven		polypropylene nonwoven and barrier laminate		cotton-like nonwoven and barrier laminate	
colour versions								
sizes	M-XXL							
additional versions	with towels and additionally wrapped in crepe paper							

\* PLUS line - elongated gowns, resistant to agents containing alcohol

Medidress

REUSABLE

compliance with EN 13795 standard

safe and functional design


high barrier to liquid penetration

high absorption capacity  
in the critical product area

high strength parameters



Reusable surgical gowns - Medidress

	Reusable surgical gowns - Medidress	
	Standard performance	High performance
name	Standard 	Perfect 
raw material colour versions	polyester fabric 	polyester fabric and barrier laminate 
sizes	S-XXL	

## Matodress



### DISPOSABLE

compliance with EN 13795 standard

antistatic properties


high strength parameters

water vapour permeability

attractive and practical design

wide range of sizes

#### Disposable surgical clothing - Matodress

	Surgical clothing set	Surgical dress
name		
raw material colour versions	polypropylene nonwoven	
sizes	XS-XXXL	S-XL

Medidress

REUSABLE

compliance with EN 13795 standard

antistatic properties

high resistance


water vapour permeability

attractive and practical design

wide range of sizes



Reusable surgical clothing - Medidress

	Surgical clothing set		Surgical dress
	male	female	
name			
raw material colour versions	polyester fabric 		cotton and polyester fabric 
sizes	XS-3XL / 158-164; 170-176; 182-188		



 **TZMO SA**  
| [tzmo-global.com](http://tzmo-global.com)  
| [matopat-global.com](http://matopat-global.com)