



Reusable bed linen

Bed linen, as a basic product used in hospitals, hotels and other institutions, must be of very high quality due to constant turnover.

They say "from thread to fabric", and this is precisely how the bed linen we offer is produced.

We produce fabrics starting from the yarn itself, so from the very beginning, carefully selected raw materials are under the watchful eye of quality control.

We make sure that our bed linen fabrics not only have

very good strength properties



but also look good and provide comfort to patients.

The fabrics from which the bed linen is made **are certified in accordance with textile standards**

for the healthcare system and hold the OEKO-TEX® STANDARD 100 certification.

Bed linen

with a durable
bacteriostatic finish




- inhibits the activity and multiplication of Staphylococcus aureus bacteria
- eliminates unpleasant sweat odours
- has anti-allergic properties

Medilinen is a line of products that perfectly meets the quality and utility requirements. It is intended for use during the hospitalisation of patients in medical facilities as a product designed to protect the patient against infections.

Medilinen bed linen meets the requirements for a medical device in accordance with Regulation (EU) 2017/745 of the European Parliament and of the Council (MDR).


Pillow case with a flap

- 1 length - 40 cm
- 2 width - 40 cm
- 3 flap width - 20 cm
- 4 length from edge to fastening - 5 cm
- 5 opening length - 30 cm




Pillow case with a flap

- 1 length - 70 cm
- 2 width - 80 cm
- 3 flap width - 20 cm
- 4 length from edge to fastening - 5 cm
- 5 opening length - 65 cm




Duvet cover with a flap

- 1 length - 210 cm
- 2 width - 160 cm
- 3 flap width - 30 cm
- 4 length from edge to fastening - 5 cm
- 5 opening length - 145-150 cm




Underpad

- 1 length - 160 cm
- 2 width - 90 cm




Sheet

- 1 length - 250 cm
- 2 width - 160 cm




Double-sided pillow case

- 1 length - 70 cm
- 2 width - 80 cm
- 3 length from edge to fastening - 10 cm




Double-sided duvet cover

- 1 length - 210 cm
- 2 width - 160 cm
- 3 length from edge to fastening - 20 cm



Bed linen coding




Reusable bed linen

4


- 1 Bed linen symbol:** PD - duvet cover, PM - pillow case, PR - sheet, PRG - sheet with elastic band, PK - underpad
- 2 Product size** – from the table located at the end of the catalogue
- 3 Fabric type and colour** – fabric name-colour number – from the palette on pages 5 and 6, e.g. MEDYK3-3
- 4 Information about double-sidedness** – DW (double-sided)
- 5 Information about bacteriostatic finishing** – B (bacteriostatic)

Colour palette



Attention:

The colour of the inscription woven into the fabric's structure can be any



 Fabrics for bed linen available in the CitoNet service



▲ Attention:

Makao fabric available without the inscription woven into the fabric's structure





PM/80X70/MAKAO-7

PR/280X160/MEDYK4-1

PD/210X160/MAKAO-6

Reusable bed linen





▲ Attention:

Tiwoli fabric available without the inscription woven into the fabric's structure

7

Reusable bed linen



PM/80X70/MAKAO-1/B



PM/70X50/MAKAO-1/B


PR/280X160/MEDYK3-1/B



PD/210X160/MAKAO-3/B


Reusable bed linen





▲ Attention:

Medyk3 fabric with decorative stripes available only for CitoNet - service



Double-sided bed linen

PM/60X50/MEDYK3-2/B




PRG/220X90/MEDYK3-5/DW

PK/150X90/MEDYK-1/B



PD/200X140/MEDYK3-2/B





▲ Attention:

The minimum order for fabric with a coloured interwoven inscription is 100 sets of bed linen


Reusable bed linen



PM/80X70/MEDYK3-3/DW



PR/280X160/MEDYK-5/DW



PK/160X90/MEDYK-3

PD/210X160/MEDYK3-3/DW

Reusable bed linen



Double-sided bed linen

PM/70X50/NORIS-4



PD/140X115/NORIS-4

PGR/220X90/MEDYK3-2

15

Reusable bed linen



Positive Opinion
of the Institute
of Mother and Child
no. OP-5321-5322.

Opinions can be checked on the website
www.epozytywnaopinia.pl

PM/62X50/NORIS-1



PR/250X160/MEDYK-2

PD/140X100/NORIS-1

Reusable bed linen

Positive Opinion
of the Institute
of Mother and Child
no. OP-5321-5322.

Opinions can be checked on the website
www.epozytywnaopinia.pl



MEDYK

1

FABRIC / BASIC WEIGHT	MEDYK / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	BLUE
FINISH	fabric dyed with VAT dyes

2

FABRIC / BASIC WEIGHT	MEDYK / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	WHITE
FINISH	bleached fabric

3

FABRIC / BASIC WEIGHT	MEDYK / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	GREEN
FINISH	fabric dyed with VAT dyes

MEDYK3

1

FABRIC / BASIC WEIGHT	MEDYK3 / 140g/m ²
RAW MATERIAL	50% COTTON 50% POLYESTER
COLOUR	WHITE
FINISH	bleached fabric

2

FABRIC / BASIC WEIGHT	MEDYK3 / 140g/m ²
RAW MATERIAL	50% COTTON 50% POLYESTER
COLOUR	BLUE
FINISH	fabric dyed with VAT dyes

3

FABRIC / BASIC WEIGHT	MEDYK3 / 140g/m ²
RAW MATERIAL	50% COTTON 50% POLYESTER
COLOUR	WHITE with a decorative graphite stripe
FINISH	bleached fabric

4

FABRIC / BASIC WEIGHT	MEDYK3 / 140g/m ²
RAW MATERIAL	50% COTTON 50% POLYESTER
COLOUR	WHITE with coloured inscription
FINISH	bleached fabric

5

FABRIC / BASIC WEIGHT	MEDYK3 / 140g/m ²
RAW MATERIAL	50% COTTON 50% POLYESTER
COLOUR	WHITE with a decorative grey stripe
FINISH	bleached fabric

6

FABRIC / BASIC WEIGHT	MEDYK3 / 140g/m ²
RAW MATERIAL	50% COTTON 50% POLYESTER
COLOUR	BLUE with coloured inscription
FINISH	fabric dyed with VAT dyes

MEDYK4

1

FABRIC / BASIC WEIGHT	MEDYK4 / 120g/m ²
RAW MATERIAL	65% COTTON 35% POLYESTER
COLOUR	WHITE with a graphite inscription
FINISH	bleached fabric

TIWOLI

1

FABRIC / BASIC WEIGHT	TIWOLI / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	LIGHT BLUE
FINISH	fabric dyed with VAT dyes

2


FABRIC / BASIC WEIGHT	TIWOLI / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	WHITE
FINISH	bleached fabric

3

FABRIC / BASIC WEIGHT	TIWOLI / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	ECRU
FINISH	fabric dyed with VAT dyes


MAKAO

1	
FABRIC / BASIC WEIGHT	MAKAO / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	BURGUNDY
FINISH	fabric dyed with VAT dyes

2	
FABRIC / BASIC WEIGHT	MAKAO / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	ECRU
FINISH	fabric dyed with VAT dyes

3	
FABRIC / BASIC WEIGHT	MAKAO / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	LIGHT BLUE
FINISH	fabric dyed with VAT dyes

4	
FABRIC / BASIC WEIGHT	MAKAO / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	GRAY
FINISH	fabric dyed with VAT dyes

5	
FABRIC / BASIC WEIGHT	MAKAO / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	WHITE
FINISH	fabric dyed with VAT dyes

6	
FABRIC / BASIC WEIGHT	MAKAO / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	LIGHT GREY
FINISH	fabric dyed with VAT dyes

7	
FABRIC / BASIC WEIGHT	MAKAO / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	LILAC
FINISH	fabric dyed with VAT dyes

8	
FABRIC / BASIC WEIGHT	MAKAO / 160g/m ²
RAW MATERIAL	100% COTTON
COLOUR	BEIGE
FINISH	fabric dyed with VAT dyes

NORIS

1	
FABRIC / BASIC WEIGHT	NORIS / 150g/m ²
RAW MATERIAL	100% COTTON
COLOUR	CHILDREN'S PATTERNS
FINISH	fabric with printed children's patterns

2	
FABRIC / BASIC WEIGHT	NORIS / 150g/m ²
RAW MATERIAL	100% COTTON
COLOUR	CHILDREN'S PATTERNS
FINISH	fabric with printed children's patterns

3	
FABRIC / BASIC WEIGHT	NORIS / 150g/m ²
RAW MATERIAL	100% COTTON
COLOUR	CHILDREN'S PATTERNS with coloured inscription
FINISH	fabric with printed children's patterns

4	
FABRIC / BASIC WEIGHT	NORIS / 150g/m ²
RAW MATERIAL	100% COTTON
COLOUR	CHILDREN'S PATTERNS with coloured inscription
FINISH	fabric with printed children's patterns



Standard bed linen sizes and possibility of interweaving the inscription

	Medyk	Medyk 3	Medyk 4	Tiwoli	Makao	Noris
Duvet cover - PD						
210 x 160	⊙	⊙	⊙	⊙	⊙	⊙
210 x 160 double-sided		⊙				
200 x 140	⊙	⊙	⊙	⊙	⊙	
140 x 115						⊙
140 x 100						⊙
160 x 100						⊙
85 x 75						⊙
Pillow case - PM						
80 x 70	⊙	⊙	⊙	⊙	⊙	⊙
80 x 70 double-sided		⊙				
70 x 50				⊙	⊙	⊙
62 x 50						⊙
60 x 50						⊙
40 x 40			⊙	⊙	⊙	
Underpad - PK						
160 x 90	⊙	⊙	⊙			
160 x 90 double-sided		⊙				
150 x 90	⊙	⊙	⊙			
Sheet - PR						
280 x 160	⊙	⊙	⊙	⊙	⊙	⊙
280 x 160 double-sided		⊙				
260 x 160	⊙	⊙	⊙	⊙	⊙	⊙
250 x 160	⊙	⊙	⊙	⊙	⊙	⊙
240 x 160	⊙	⊙	⊙	⊙	⊙	⊙
210 x 160	⊙	⊙	⊙	⊙	⊙	⊙
220 x 90 with elastic band	⊙	⊙	⊙	⊙	⊙	⊙
160 x 100	⊙	⊙	⊙			⊙
145 x 120						⊙
120 x 70						⊙

Dimensions given in [cm]

- ⊙ no possibility of interweaving the inscription
- ⊙ possibility of interweaving the inscription
- ⊙ with a decorative graphite stripe
- ⊙ with a decorative gray stripe
- fabric for CitoNet