
Medisorb SILVER PAD Silver dressing

Managing the infection

- the absorbent pad contains activated metallic silver which has antibacterial properties
- the dressing absorbs excess secretion and maintains a moist wound environment
- the elastic nonwoven fabric allows the product to adhere to the contours of the body shape
- steam permeable, breathable – allows for gas exchange



Application:

- hard-to-heal wounds
- infected wounds and clean wounds as a dressing to prevent infections
- graft donor sites

Types:

- sterile

Sizes and packaging:

- 10 cm x 10 cm, packed by 5 pcs
- 10 cm x 20 cm, packed by 5 pcs



MD

Sterile devices are certified by the Notified Body Superiore di Sanità and bear the CE mark with number 0373

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026