Matodrape/Camera Cover Film covers for medical device cables

Made of transparent polyethylene film.
The elasticity and softness of the material allows for easy and safe manipulation during surgery.
It is equipped with an adhesive tape and a cardboard box with cut-outs, so that the cover can be easily fixed on the device.

Application:

- to cover medical devices, universal for all cameras and endoscopes, fibre optic and cables

Size and packing:

- sizes: 13 cm x 235 cm, 13 cm x 244 cm, 15 cm x 250 cm
- sterile (sterilised with radiation or ethylene oxide): packed in a paper-film by 1 pc





MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026