Ambulex & Ambulex P Examination and disposable protective gloves

Combination of comfort and safety.
Latex examination gloves are characterised by high elasticity and mechanical strength.
Provide an effective protective barrier, and a beaded cuff enables efficient and secure donning.

Ambulex – examination and disposable protective, powdered

- powdered with natural corn starch, which insolates the hand from latex layer
- textured fingertips provides a better grip for efficient instrument handling
- ambidextrous – fit both hands

Ambulex gloves classified as personal protective equipment (PPE) of category III, comply with EN ISO 374. Before using the gloves refer to the add. information for users. For detailed information on chemical resistance, please contact your medical sales representative.



Ambulex P – examination and disposable protective gloves, powder-free – especially recommended for people with sensitive skin

- made from latex, powder-free
- polymer coating protects enables effortless damp and dry hand donning
- rough surface provides a secure grip for efficient instrument handling in dry and wet conditions
- ambidextrous – fit both hands

Ambulex P gloves classified as personal protective equipment (PPE) of category III, comply with EN ISO 374. Before using the gloves refer to the add. information for users. For detailed information on chemical resistance, please contact your medical sales representative.

Ambulex P – Declaration of Conformity

Application:

- widespread use in healthcare
- reduce the risk of infection during the study and treatment, and contact with contaminated medical materials

Types:

- non sterile
- powdered and powder-free

Size and packing:

- XS, S, M, L, XL in a dispenser box by 100 pcs
- powdered Ambulex gloves size M packed in a box by 10 pcs dedicated for household use





MD

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026