Matodrape Medical fixing tape

The elasticity of the raw material, makes it easy to adjust to any shape, and the applied “fingerlifts” allow for quick application.
The Velcro tape is composed of two parts: hooks and loops. Both parts have a glue layer and a protective release paper with “fingerlifts”. One part of the tape is to be attached to a given place on a drape, and the second part can be removed or attached repeatedly again as often as needed.

Application:

- for fixing drapes, tubes, cables


Types:

- with or without Velcro

Size and packing:

- sterile (sterilised with ethylene oxide): packed in a paper-film pouch by 1 or 2 pcs
- sizes: 50 cm x 9 cm, 50 cm x 10 cm, 90 cm x 5 cm, 75 cm x 10 cm (without Velcro) and 25 cm x 2 cm (with Velcro)





MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026