Medisorb A Calcium alginate dressing

Medisorb A is a sterile, soft, and highly absorbent dressing made from 100% calcium alginate fibres.

Effect

Properties of calcium alginate dressing:
– calcium alginate fibres, which the dressing is made from, react with wound exudate to form a soft, moist gel, which maintains a moist wound healing environment
– the moist gel does not stick to the wound, reducing the risk of injury during dressing change
– for deep wounds, the dressing can easily adapt to the shape of the wound
– calcium alginate fibres release calcium into the wound, which helps to stop minor bleeding
– the integrity of the wet dressing allows it to be removed in one piece during dressing change

Indications

Medisorb A is intended for the treatment of superficial or deep wounds with moderate to heavy exudate, such as pressure sores, leg ulcers, cavitary wounds, and skin graft sites.
The dressing can stay in place for up to 7 days, depending on the patient's condition and the amount of exudate.

Directions for use

Clean and disinfect the wound before applying the dressing. It requires secondary dressing protection.

Additional remarks

Latex-free dressing.
Medisorb A is intended for professional use only.

Contraindications

Contraindications: Do not use for surgical implants, third-degree burns, or dry necrotic wounds.

Available sizes


10 cm x 10 cm – 10 pieces

Distributor

Toruńskie Zakłady Materiałów Opatrunkowych SA

Manufacturer

Eurofarm S.p.A., Italy





MD

Sterile devices are certified by the Notified Body Superiore di Sanità and bear the CE mark with number 0373

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026