Ambulex NITRYL Nitrile examination gloves

Thinking of people who are allergic to latex.
Nitrile examination gloves are characterised by the absence of natural rubber latex and high tensile strength.
Provide an effective protective barrier, and a beaded cuff enables efficient and secure donning.



Ambulex NITRYL – nitrile examination gloves, powder-free

- powder-free
- high tensile strength both at the beginning and end of the validity period
- textured fingertips provide a secure grip for efficient instrument handling in dry and wet conditions
- ambidextrous – fit both hands
- available colours: white, violet, black (especially for tattoo studios), elongated white (an elongated cuff allows for better protection against infections)

The powder-free nitrile gloves are classified as personal protective equipment of category III, comply with EN ISO 374. Before using the gloves refer to the add. information for users. For detailed information on chemical resistance, please contact your medical sales representative.

Application:

- widespread use in healthcare
- reduce the risk of infection during examinations and procedures, and contact with contaminated medical materials





Types:

- non sterile
- powder-free

Size and packing:

- white, violet, elongated white: XS, S, M, L, XL
- black: S, M, L
- in a dispenser box by 100 pcs

MD

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026