Matodrape Medical equipment cover

Made of tear resistant SMS nonwoven. The covers are fixed on a device by means of ties or elastic bands.

Application:


- to cover medical devices


Types:


- one-piece C-arm cover (A, B)
- two-piece C-arm cover (C)

Size and packing:

- sterile (steam sterilised): packed in a paper-film pouch by 1 pc







MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026