Fixopore Nonwoven plaster with absorbent pad

- made of white elastic nonwoven, which easily adjusts to the body shape
- porous fabric is skin-friendly and breathable
- centrally placed absorbent pad is secured with micro-net to protect the dressing from sticking to the wound
- acrylic adhesive does not cause skin irritation
- the length of the plaster enables to adjust the size to individual needs

Application:

- for dressing abrasions, minor cuts and wounds in people with sensitive skin

Types:

- non sterile

Size and packing:

- available sizes: 6 cm x 1 m, 8 cm x 1 m
- non sterile: packed by 1 pc







MD

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026