Silicone dressings for scar correction Codosil™ ADHESIVE

Aesthetic scar correction.

- used for the rehabilitation of hypertrophic scars and keloids
- regularly caused reduces the size and redness of scars
- can be cut to the size and shape of the scar – so it acts only on the scar and does not cause maceration of healthy skin; the rest of the dressing may be used for further rehabilitation
- reusable – since the cut piece of the dressing can be washed it is very economical
- used in combination therapy with the Codopress® compression garment

Application:

- in the rehabilitation of scars, including overgrown scars and keloids
- for prophylactic purposes in people with a tendency to develop hypertrophic scars and keloids

Types:

- non sterile
- square and rectangular shape
- 2 shapes for breasts

Size and packing:

- non sterile: 7 cm x 7 cm, 7 cm x 14 cm, 14 cm x 14 cm, packed by 1 pc
- non sterile: Breast 1 and Breast 2, packed by 1 pc







MD

Sterile devices are certified by the Notified Body PCBC SA and bear the CE mark with number 1434

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026