- made of durable fabric that adjusts perfectly to the body
- serrated edges on both sides for easier tearing off the desired length in the hand
- the adhesive contains zinc oxide, which lowers the risk of skin irritation

Application:


- for supporting dressings
- for fixing probes, catheters and heavier medical devices

Types:

- non sterile

Size and packing:

- available length: 5 m and a width of 1.25 cm, 2.5 cm and 5 cm
- non sterile: packaging depending on the size








MD

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026