Matolast UNIVERSAL Elastic bandage

A delicate elastic bandage with a high stretch of 150% and universal use.

Reusable, retains constant elasticity, even after several washings.

Each bandage comes with a clip that facilitates dressing attachment.

Application:

- in the prevention and treatment of musculoskeletal injuries, such as sprains or strains
- as a light compression bandage (for the treatment of superficial venous system)
- as a support bandage – for fixing dressings

Types:

- non sterile

Colours:

- red
- green
- beige

Size and packing:

- non sterile 8 cm, 10 cm, 12 cm, 15 cm length 4 m film packaging 1 szt. – Length measured in a stretched condition.

Length measured in a stretched condition.







MD

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026