Matovis Support bandage

Light support non-elastic bandage designed for a variety of applications.
High comfort of use is provided by high breathability and evenly finished and non-shredding edges.

Application:

- as a supporting bandage – for fixing dressings
- for fixing and supporting medical devices

Types:

- sterile (steam sterilised) and non sterile

Size and packing:

- sterile: 5 cm, 10 cm and 15 cm and a length of 4 m, packed by 1 pc
- sterile: 10 cm and 15 cm and a length of 4 m, packed by 2 pc
- non sterile: 5 cm, 10 cm, 15 cm and a length of 4 m, packed in film packaging by 1 pc







MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>

©Copyright 2026 en.matopat-global.com | Last updated: 04.11.2026